



# BETAJECT® Lite

## How to get your BETAJECT® Lite

To order the BETAJECT® Lite, please contact your local Bayer Schering Pharma AG organisation or Betaferon® nurse.

For further information on Betaferon®, please refer to our MS Information Gateway website at [www.ms-gateway.com](http://www.ms-gateway.com).



EU2006.0816



The new generation of autoinjectors

*Suitable for the new Betaferon® injection system. Autoinjectors have been proven to reduce injection site reactions. They are particularly suitable for people with a fear of needles.*



## The new BETAJECT® Lite

### Customised

- Improved design with new ergonomic grip
- Constant spring force throughout the whole injection process for continuous release
- Comes in a nice case

### Flexible

- Flexible depth adjustment possible  
Please ask your health care professional for advice
- Suitable for dose deliveries from 0.25 to 1.0 ml
- Can be used with the thinnest needle size available (30 G)

### Convenient

- Insert the syringe with the needle cap still attached
- 3 large optical markers indicate the end of injection
- Prolonged injection time for gentle administration

## Ready, Steady, Go

### Step 1: Prepare the device



### Step 2: Insert the ready-to-use syringe



### Step 3: Inject



The smart autoinjector with advanced technology.